**Patient:** HAKAM SINGH 53Y/M (male, \*1968-12-11, #59262)

Referring Physician: DR ANURAG LAMBA
Study: CEMRI BRAIN (#1)
Series: PhoenixZIPReport (#99)

Protocol: Phoenix Document

Manufacturer: SIEMENS (Spectra)
Completion Flag: PARTIAL

Verification Flag: UNVERIFIED

**Content Date/Time:** 2021-12-11 14:39:34

## **Phoenix ZIP Acquisition Report**

Concept Modifier: Document Title Modifier = Phoenix ZIP Document (MR-20, 99SMS\_CTMR)

Concept Modifier: Language of Content Item and Descendents = English (eng, ISO639\_2)

Observation Context: Observer Type = Person (121006, DCM)

Observation Context: Observer Name = n.a.

## **Current Procedure Description**

**Modality:** 

Magnetic Resonanz (MR, DCM)

**Performed Procedure Step:** 

n.a.

## **Acquisition Protocol**

## **Protocol Header**

**System Type:** 

Spectra

**Software Version:** 

syngo MR E11

Field Strength:

2.89362

**Protocol Name:** 

t2\_spc\_tra\_iso

**Sequence Name:** 

tse\_vfl

**Acquisition Time:** 

00:04:34

**Orientation:** 

158

**Resolution Matrix:** 

384x384

**Number of Slices:** 

1 1

FoV Read:

200 mm

FoV Phase:
100 %
Slice Thickness:
0.95 mm
Slice Gap:
empty [Measurement not attempted]
Repetition Time:
1000 msec
Echo Time:
138 msec
Flip Angle:
120 deg
Fat Suppression:
35
Flow Compensation:
301
Water Suppression:
Subtraction:
False St. 114 1
First Signal Mode:
236
Contrast Agent:
n.a.
Image Acquired
Acquisition Context: Image Type = "n.a."
Comment:
n.a.

This page was generated from a DICOM Structured Reporting document by OFFIS DCMTK 3.6.6.