Patient: HAKAM SINGH 53Y/M (male, *1968-12-11, #59262)

Referring Physician: DR ANURAG LAMBA
Study: CEMRI BRAIN (#1)
Series: PhoenixZIPReport (#99)

Protocol: Phoenix Document **Manufacturer:** SIEMENS (Spectra)

Completion Flag: PARTIAL
Verification Flag: UNVERIFIED

Content Date/Time: 2021-12-11 14:32:24

Phoenix ZIP Acquisition Report

Concept Modifier: Document Title Modifier = Phoenix ZIP Document (MR-20, 99SMS_CTMR)

Concept Modifier: Language of Content Item and Descendents = English (eng, ISO639_2)

Observation Context: Observer Type = Person (121006, DCM)

Observation Context: Observer Name = n.a.

Current Procedure Description

Modality:

Magnetic Resonanz (MR, DCM)

Performed Procedure Step:

n.a.

Acquisition Protocol

Protocol Header

System Type:

Spectra

Software Version:

syngo MR E11

Field Strength:

2.89362

Protocol Name:

t2_tirm_tra_dark-fluid

Sequence Name:

tse

Acquisition Time:

00:02:26

Orientation:

158

Resolution Matrix:

320x203

Number of Slices:

24 1

FoV Read:

240 mm

FoV Phase:	
90.6 %	
Slice Thickness:	
5 mm	
Slice Gap:	
30 %	
Repetition Time:	
6000 msec	
Echo Time:	
77 msec	
Flip Angle:	
150 deg	
Fat Suppression:	
33,615	
Flow Compensation:	
301	
Water Suppression:	
39	
Subtraction:	
False	
First Signal Mode:	
236	
Contrast Agent:	
n.a.	
Image Acquired	
•	
Acquisition Context: Image Type = "n.a."	
Comment:	
n.a.	

This page was generated from a DICOM Structured Reporting document by OFFIS DCMTK 3.6.6.