**Patient:** HAKAM SINGH 53Y/M (male, \*1968-12-11, #59262)

Referring Physician: DR ANURAG LAMBA
Study: CEMRI BRAIN (#1)
Series: PhoenixZIPReport (#99)

**Protocol:** Phoenix Document **Manufacturer:** SIEMENS (Spectra)

Completion Flag: PARTIAL

**Verification Flag:** 

**Content Date/Time:** 2021-12-11 14:35:47

## **Phoenix ZIP Acquisition Report**

**UNVERIFIED** 

Concept Modifier: Document Title Modifier = Phoenix ZIP Document (MR-20, 99SMS\_CTMR)

Concept Modifier: Language of Content Item and Descendents = English (eng, ISO639\_2)

Observation Context: Observer Type = Person (121006, DCM)

Observation Context: Observer Name = n.a.

## **Current Procedure Description**

**Modality:** 

Magnetic Resonanz (MR, DCM)

**Performed Procedure Step:** 

n.a.

## **Acquisition Protocol**

## **Protocol Header**

**System Type:** 

Spectra

**Software Version:** 

syngo MR E11

**Field Strength:** 

2.89362

**Protocol Name:** 

t2\_tse\_SAG

**Sequence Name:** 

tse

**Acquisition Time:** 

00:00:47

**Orientation:** 

150

**Resolution Matrix:** 

384x278

**Number of Slices:** 

20 1

**FoV Read:** 

260 mm

FoV Phase:
90.6 %
Slice Thickness:
4.5 mm
Slice Gap:
50 %
Repetition Time:
4500 msec
Echo Time:
78 msec
Flip Angle:
150 deg
Fat Suppression:
35
Flow Compensation:
301
Water Suppression:
39
Subtraction:
False
First Signal Mode:
236
Contrast Agent:
n.a.
Image Acquired
Acquisition Context: Image Type = "n.a."
Comment:
n.a.

This page was generated from a DICOM Structured Reporting document by OFFIS DCMTK 3.6.6.